

For Research Use Only

RAB37 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83740-1-PBS

Basic Information

Catalog Number:	83740-1-PBS	GenBank Accession Number:	BC016615	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	326624	CloneNo.:	240838A6
Source:	Rabbit	UNIPROT ID:	Q96AX2		
Isotype:	IgG	Full Name:	RAB37, member RAS oncogene family		
Immunogen Catalog Number:	AG3717	Calculated MW:	223 aa, 25 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00690-2, RAB37 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83740-1-PBS. Detection antibody: 83740-2-PBS. Standard: Ag3717. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00690-3, RAB37 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83740-3-PBS. Detection antibody: 83740-1-PBS. Standard: Ag3717. Range: 0.625-80 ng/mL