

For Research Use Only

Rat C-reactive protein Recombinant antibody, PBS Only (Capture)

Catalog Number: 83741-1-PBS



Basic Information

Catalog Number:

83741-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_017096.3

GeneID (NCBI):

25419

UNIPROT ID:

P48199

Full Name:

C-reactive protein, pentraxin-related

Calculated MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

240759D4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00692-1, RAT C-Reactive Protein/CRP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83741-1-PBS. Detection antibody: 83741-2-PBS. Standard: Eg0942. Range: 0.39-50 ng/mL