

For Research Use Only

# IL-6 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83747-4-PBS



## Basic Information

Catalog Number:

83747-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC015511

GeneID (NCBI):

3569

ENSEMBL Gene ID:

ENSG00000136244

UNIPROT ID:

P05231

Full Name:

interleukin 6 (interferon, beta 2)

Calculated MW:

212 aa, 24 kDa

Purification Method:

Protein A purification

CloneNo.:

240670E1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00687-1, IL-6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83747-2-PBS. Detection antibody: 83747-4-PBS. Standard: Eg0117. Range: 0.625-80 ng/mL