

For Research Use Only

CKMT1A Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83750-1-PBS

Basic Information

Catalog Number: 83750-1-PBS	GenBank Accession Number: BC001926	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 548596	CloneNo.: 240825A2
Source: Rabbit	UNIPROT ID: P12532	
Isotype: IgG	Full Name: creatine kinase, mitochondrial 1A	
Immunogen Catalog Number: AG7583	Calculated MW: 47 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00706-1, CKMT1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83750-2-PBS. Detection antibody: 83750-1-PBS. Standard: Ag7583. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00706-2, CKMT1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83750-2-PBS. Detection antibody: 83750-3-PBS. Standard: Ag7583. Range: 0.313-40 ng/mL