

For Research Use Only

OPN1SW Recombinant antibody, PBS Only (Capture)



Catalog Number: 83754-2-PBS

Basic Information

Catalog Number:

83754-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20445

GenBank Accession Number:

BC156719

GeneID (NCBI):

611

UNIPROT ID:

P03999

Full Name:

opsin 1 (cone pigments), short-wave-sensitive

Calculated MW:

348 aa, 39 kDa

Purification Method:

Protein A purification

CloneNo.:

240714H5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00721-2, OPN1SW Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83754-2-PBS. Detection antibody: 83754-3-PBS. Standard: Ag20445. Range: 0.78-100 ng/mL.