

For Research Use Only

# BCL9 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83760-4-PBS



## Basic Information

Catalog Number:

83760-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19158

GenBank Accession Number:

BC116451

GeneID (NCBI):

607

UNIPROT ID:

O00512

Full Name:

B-cell CLL/lymphoma 9

Calculated MW:

1426 aa, 149 kDa

Purification Method:

Protein A purification

CloneNo.:

240682E8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926


E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)


This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.



## Selected Validation Data



Cytometric bead array standard curve of MP00739-3, BCL9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83760-3-PBS. Detection antibody: 83760-4-PBS. Standard: Ag19158. Range: 0.625-40 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83760-4-PBS against Human BCL9 was performed. The affinity constant is 0.816 nM.