

For Research Use Only

PAK4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83761-1-PBS

Basic Information

Catalog Number:	83761-1-PBS	GenBank Accession Number:	BC025282	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	10298	CloneNo.:	240686A9
Source:	Rabbit	UNIPROT ID:	O96013		
Isotype:	IgG	Full Name:	p21 protein (Cdc42/Rac)-activated kinase 4		
Immunogen Catalog Number:	AG5981	Calculated MW:	64 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00743-2, PAK4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83761-4-PBS. Detection antibody: 83761-1-PBS. Standard: Ag5981. Range: 2.5-40 ng/mL.