

For Research Use Only

# PAK4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83761-2-PBS

## Basic Information

Catalog Number:	83761-2-PBS	GenBank Accession Number:	BC025282	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	10298	CloneNo.:	240686E1
Source:	Rabbit	UNIPROT ID:	O96013		
Isotype:	IgG	Full Name:	p21 protein (Cdc42/Rac)-activated kinase 4		
Immunogen Catalog Number:	AG5981	Calculated MW:	64 kDa		

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00743-1, PAK4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83761-3-PBS. Detection antibody: 83761-2-PBS. Standard: Ag5981. Range: 1.25-80 ng/mL.