

For Research Use Only

PAK4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83761-4-PBS

Basic Information

Catalog Number: 83761-4-PBS	GenBank Accession Number: BC025282	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10298	CloneNo.: 240686F11
Source: Rabbit	UNIPROT ID: O96013	
Isotype: IgG	Full Name: p21 protein (Cdc42/Rac)-activated kinase 4	
Immunogen Catalog Number: AG5981	Calculated MW: 64 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00743-2, PAK4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83761-4-PBS. Detection antibody: 83761-1-PBS. Standard: Ag5981. Range: 2.5-40 ng/mL.