

For Research Use Only

# EPO/Erythropoietin Recombinant antibody, PBS Only (Detector)

Catalog Number: 83762-2-PBS



## Basic Information

Catalog Number:

83762-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC093628

GeneID (NCBI):

2056

UNIPROT ID:

P01588

Full Name:

erythropoietin

Calculated MW:

193 aa, 21 kDa

Purification Method:

Protein A purification

CloneNo.:

240773D3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00751-1, EPO/Erythropoietin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83762-3-PBS. Detection antibody: 83762-2-PBS. Standard: Eg0197. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00751-2, EPO/Erythropoietin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83762-5-PBS. Detection antibody: 83762-2-PBS. Standard: Eg0197. Range: 0.313-40 ng/mL