

For Research Use Only

# HIF-1 alpha Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 83764-1-PBS



## Basic Information

Catalog Number:

83764-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0939

GenBank Accession Number:

NM\_001530.4

GeneID (NCBI):

3091

UNIPROT ID:

Q16665-1

Full Name:

hypoxia inducible factor 1, alpha subunit (basic helix-loop-helix transcription factor)

Calculated MW:

93 kDa

Purification Method:

Protein A purification

CloneNo.:

240767E7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00762-1, HIF-1 alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83764-2-PBS. Detection antibody: 83764-1-PBS. Standard: Eg0939. Range: 0.625-80 ng/mL.