

For Research Use Only

HNRNPA2B1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83773-5-PBS



Basic Information

Catalog Number:

83773-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6563

GenBank Accession Number:

BC000506

GeneID (NCBI):

3181

UNIPROT ID:

P22626

Full Name:

heterogeneous nuclear
ribonucleoprotein A2/B1

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

240831G11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00768-1, HNRNPA2B1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83773-4-PBS. Detection antibody: 83773-5-PBS. Standard: Ag6563. Range: 0.313-40 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83773-5-PBS against Human HNRNPA2B1 was performed. The affinity constant is below 1 pM.