

For Research Use Only

HNRNPA2B1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83773-8-PBS



Basic Information

Catalog Number:

83773-8-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6563

GenBank Accession Number:

BC000506

GeneID (NCBI):

3181

UNIPROT ID:

P22626

Full Name:

heterogeneous nuclear
ribonucleoprotein A2/B1

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

240831B3

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00768-4, Human HNRNPA2B1 Recombinant Matched Antibody Pair - PBS only. 83773-8-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag6563. 83773-2-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL



Serum of sixteen individual healthy human donors was measured. The HNRNPA2B1 concentration of detected samples was determined to be 1.71 ng/mL with a range of 0.16 - 5.49 ng/mL