

For Research Use Only

# ASB2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83775-1-PBS

## Basic Information

Catalog Number: 83775-1-PBS	GenBank Accession Number: BC032354	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51676	CloneNo.: 240722C8
Source: Rabbit	UNIPROT ID: Q96Q27	
Isotype: IgG	Full Name: ankyrin repeat and SOCS box- containing 2	
Immunogen Catalog Number: AG6270	Calculated MW: 65 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00730-2, ASB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83775-1-PBS. Detection antibody: 83775-2-PBS. Standard: Ag6270. Range: 3.13-100 ng/mL