

For Research Use Only

ASB2 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83775-2-PBS

Basic Information

Catalog Number:

83775-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6270

GenBank Accession Number:

BC032354

GeneID (NCBI):

51676

UNIPROT ID:

Q96Q27

Full Name:

ankyrin repeat and SOCS box-containing 2

Calculated MW:

65 kDa

Purification Method:

Protein A purification

CloneNo.:

240722G1

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00730-1, ASB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83775-2-PBS. Detection antibody: 83775-3-PBS. Standard: Ag6270. Range: 1.56-100 ng/mL



Cytometric bead array standard curve of MP00730-2, ASB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83775-1-PBS. Detection antibody: 83775-2-PBS. Standard: Ag6270. Range: 3.13-100 ng/mL