

For Research Use Only

LILRA4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83777-1-PBS



Basic Information

Catalog Number:

83777-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15721

GenBank Accession Number:

BC109198

GeneID (NCBI):

23547

UNIPROT ID:

P59901

Full Name:

leukocyte immunoglobulin-like
receptor, subfamily A (with TM
domain), member 4

Calculated MW:

499 aa, 55 kDa

Purification Method:

Protein A purification

CloneNo.:

240820B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00735-1, LILRA4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83777-2-PBS. Detection antibody: 83777-1-PBS. Standard: Ag15721. Range: 6.25-200 ng/mL