

For Research Use Only

CTLA-4/CD152 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83778-5-PBS



Basic Information

Catalog Number:

83778-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC070162

GeneID (NCBI):

1493

UNIPROT ID:

P16410

Full Name:

cytotoxic T-lymphocyte-associated
protein 4

Calculated MW:

223 aa, 25 kDa

Purification Method:

Protein A purification

CloneNo.:

240793F7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00736-1, CTLA4/CD152 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83778-1-PBS. Detection antibody: 83778-5-PBS. Standard: Eg0972. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00736-2, CTLA4/CD152 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83778-4-PBS. Detection antibody: 83778-5-PBS. Standard: Eg0972. Range: 0.313-40 ng/mL