

For Research Use Only

# BAG3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83779-3-PBS



## Basic Information

Catalog Number:

83779-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0956

GenBank Accession Number:

BC006418

GeneID (NCBI):

9531

UNIPROT ID:

O95817

Full Name:

BCL2-associated athanogene 3

Calculated MW:

61 kDa

Purification Method:

Protein A purification

CloneNo.:

240860D4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00740-2, BAG3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83779-1-PBS. Detection antibody: 83779-3-PBS. Standard: Ag0956. Range: 0.625-40 ng/mL