

For Research Use Only

IL21 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83781-1-PBS



Basic Information

Catalog Number:

83781-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC066258

GeneID (NCBI):

59067

UNIPROT ID:

Q9HBE4

Full Name:

interleukin 21

Calculated MW:

162 aa, 19 kDa

Purification Method:

Protein A purification

CloneNo.:

240718A8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00745-1, IL-21 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83781-1-PBS. Detection antibody: 83781-2-PBS. Standard: Eg0633. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83781-1-PBS against Human IL21 was performed. The affinity constant is 7.06 nM.