

For Research Use Only

MAFB Recombinant antibody, PBS Only (Detector)



Catalog Number: 83783-3-PBS

Basic Information

Catalog Number:	83783-3-PBS	GenBank Accession Number:	BC036689	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	9935	CloneNo.:	240821C12
Source:	Rabbit	ENSEMBL Gene ID:	ENSG00000204103		
Isotype:	IgG	UNIPROT ID:	Q9Y5Q3		
Immunogen Catalog Number:	AG14102	Full Name:	v-maf musculoaponeurotic fibrosarcoma oncogene homolog B (avian)		
		Calculated MW:	323 aa, 36 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00750-3, MAFB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83783-2-PBS. Detection antibody: 83783-3-PBS. Standard: Ag14102. Range: 0.625-40 ng/mL