

For Research Use Only

ATOX1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83785-1-PBS

Basic Information

Catalog Number:

83785-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18460

GenBank Accession Number:

BC112248

GeneID (NCBI):

475

UNIPROT ID:

O00244

Full Name:

ATX1 antioxidant protein 1 homolog (yeast)

Calculated MW:

68 aa, 8 kDa

Purification Method:

Protein A purification

CloneNo.:

240679C6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00767-2, ATOX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83785-3-PBS. Detection antibody: 83785-1-PBS. Standard: Ag18460. Range: 0.625-80 ng/mL