

For Research Use Only

# BRUNOL5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83789-2-PBS

## Basic Information

Catalog Number: 83789-2-PBS	GenBank Accession Number: BC028101	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 60680	CloneNo.: 240693G5
Source: Rabbit	UNIPROT ID: Q8N6W0	
Isotype: IgG	Full Name: bruno-like 5, RNA binding protein (Drosophila)	
Immunogen Catalog Number: AG4416	Calculated MW: 485 aa, 52 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926


E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**



# Selected Validation Data



Cytometric bead array standard curve of MP00774-1, BRUNOL5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83789-3-PBS. Detection antibody: 83789-2-PBS. Standard: Ag4416. Range: 0.625-80 ng/mL