

For Research Use Only

Mouse IL-6 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83795-2-PBS

Basic Information

Catalog Number:

83795-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_031168

GeneID (NCBI):

16193

UNIPROT ID:

P08505

Full Name:

interleukin 6

Calculated MW:

24 kDa

Purification Method:

Protein A purification

CloneNo.:

240701E6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00756-2, MOUSE IL-6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83795-4-PBS. Detection antibody: 83795-2-PBS. Standard: Eg0973. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00756-3, MOUSE IL-6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83795-5-PBS. Detection antibody: 83795-2-PBS. Standard: Eg0973. Range: 0.625-80 ng/mL