

For Research Use Only

RXRA Recombinant antibody, PBS Only (Capture)

Catalog Number: 83796-2-PBS



Basic Information

Catalog Number:

83796-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15651

GenBank Accession Number:

BC110998

GeneID (NCBI):

6256

UNIPROT ID:

P19793

Full Name:

retinoid X receptor, alpha

Calculated MW:

462 aa, 51 kDa

Purification Method:

Protein A purification

CloneNo.:

240878C1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00759-2, RXRA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83796-2-PBS. Detection antibody: 83796-1-PBS. Standard: Ag15651. Range: 0.625-80 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83796-2-PBS against Human RXRA was performed. The affinity constant is 1.20 nM.