

For Research Use Only

Mouse Cdkn2a/p16INK4a Recombinant antibody, PBS Only (Capture)

Catalog Number: 83797-4-PBS



Basic Information

Catalog Number:

83797-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28927

GenBank Accession Number:

NM_001040654

GeneID (NCBI):

12578

UNIPROT ID:

P51480

Full Name:

cyclin-dependent kinase inhibitor 2A

Calculated MW:

18 kDa

Purification Method:

Protein A purification

CloneNo.:

240935E8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00764-3, MOUSE Cdkn2a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83797-4-PBS. Detection antibody: 83797-6-PBS. Standard: Ag28927. Range: 1.563-50 ng/mL.