

For Research Use Only

# YAF2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83798-1-PBS



## Basic Information

Catalog Number:

83798-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5023

GenBank Accession Number:

BC037777

GeneID (NCBI):

10138

UNIPROT ID:

Q8IY57

Full Name:

YY1 associated factor 2

Calculated MW:

204 aa, 23 kDa

Purification Method:

Protein A purification

CloneNo.:

240834E7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00726-3, YAF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83798-5-PBS. Detection antibody: 83798-1-PBS. Standard: Ag5023. Range: 0.625-40 ng/mL