

For Research Use Only

TBX21/T-bet Recombinant monoclonal antibody, PBS Only

Catalog Number: 83804-10-PBS



Basic Information

Catalog Number:

83804-10-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4618

GenBank Accession Number:

BC039739

GeneID (NCBI):

30009

ENSEMBL Gene ID:

ENSG00000073861

UNIPROT ID:

Q9UL17

Full Name:

T-box 21

Calculated MW:

535 aa, 58 kDa

Purification Method:

Protein A purification

CloneNo.:

240849B11

Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human

Background Information

TBX21 is a transcription factor that drives the Th1 immune response primarily through promoting expression of the IFN- γ gene. TBX21 initiates TH1 lineage development from naive TH precursor cells both by activating TH1 genetic programs and by repressing the opposing TH2 programs. It has also been shown to play an influential role in inflammatory processes such as inflammatory bowel disease, autoimmune diseases such as rheumatoid arthritis, and immune-mediated conditions like asthma.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 83804-10-RR (TBX21, T-bet antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 83804-10-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 83804-10-RR (TBX21, T-bet antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 83804-10-PBS in a different storage buffer formulation.