

For Research Use Only

TBX21/T-bet Recombinant antibody, PBS Only (Detector)

Catalog Number: 83804-3-PBS



Basic Information

Catalog Number:

83804-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4618

GenBank Accession Number:

BC039739

GeneID (NCBI):

30009

ENSEMBL Gene ID:

ENSG00000073861

UNIPROT ID:

Q9UL17

Full Name:

T-box 21

Calculated MW:

535 aa, 58 kDa

Purification Method:

Protein A purification

CloneNo.:

240849B7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00748-2, TBX21,T-bet Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83804-1-PBS. Detection antibody: 83804-3-PBS. Standard: Ag4618. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00748-3, TBX21,T-bet Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83804-2-PBS. Detection antibody: 83804-3-PBS. Standard: Ag4618. Range: 0.313-40 ng/mL