

For Research Use Only

CSNK2A1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83805-2-PBS



Basic Information

Catalog Number:

83805-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1452

GenBank Accession Number:

BC011668

GeneID (NCBI):

1457

UNIPROT ID:

P68400

Full Name:

casein kinase 2, alpha 1 polypeptide

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

240862G9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00749-1, CSNK2A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83805-3-PBS. Detection antibody: 83805-2-PBS. Standard: Ag1452. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00749-2, CSNK2A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83805-2-PBS. Detection antibody: 83805-1-PBS. Standard: Ag1452. Range: 3.13-100 ng/mL