

For Research Use Only

TP53BP1-Specific Recombinant antibody, PBS Only (Detector)

Catalog Number: 83809-3-PBS



Basic Information

Catalog Number:

83809-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_005657

GeneID (NCBI):

7158

UNIPROT ID:

Q12888

Full Name:

tumor protein p53 binding protein 1

Calculated MW:

214 kDa

Purification Method:

Protein A purification

CloneNo.:

240638G10

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01611-1, Human TP53BP1 Recombinant Matched Antibody Pair - PBS only. 83809-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY00622. 83809-3-PBS was HRP conjugated as the detection antibody. Range: 1.56-50 ng/mL



The mean TP53BP1 concentration was determined to be 9.94 ng/mL in HT-29 cell extract based on a 7.70 mg/mL extract load.