

For Research Use Only

# PSENEN Recombinant antibody, PBS Only



Catalog Number: 83818-1-PBS

## Basic Information

Catalog Number: 83818-1-PBS	GenBank Accession Number: BC009575	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55851	CloneNo.: 240596H2
Source: Rabbit	UNIPROT ID: Q9NZ42	
Isotype: IgG	Full Name: presenilin enhancer 2 homolog (C. elegans)	
Immunogen Catalog Number: AG9243	Calculated MW: 101 aa, 12 kDa	

## Applications

Tested Applications:  
ELISA, WB

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83818-1-RR (PSENEN antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83818-1-PBS in a different storage buffer formulation.