

For Research Use Only

CPNE3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83820-4-PBS

Basic Information

Catalog Number:

83820-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1634

GenBank Accession Number:

BC036242

GeneID (NCBI):

8895

UNIPROT ID:

O75131

Full Name:

copine III

Calculated MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

240871G10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00783-1, CPNE3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83820-4-PBS. Detection antibody: 83820-2-PBS. Standard: Ag1634. Range: 0.625-80 ng/mL