

For Research Use Only

EXTL2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83821-1-PBS

Basic Information

| | | |
|-------------------------------------|---|--|
| Catalog Number: 83821-1-PBS | GenBank Accession Number: BC036015 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 2135 | CloneNo.: 240919B10 |
| Source: Rabbit | UNIPROT ID: Q9UBQ6 | |
| Isotype: IgG | Full Name: exostoses (multiple)-like 2 | |
| Immunogen Catalog Number: AG3548 | Calculated MW: 330 aa, 37 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00784-1, EXTL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83821-2-PBS. Detection antibody: 83821-1-PBS. Standard: Ag3548. Range: 0.313-20 ng/mL