

For Research Use Only

# ZFP64 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83825-2-PBS



## Basic Information

Catalog Number:

83825-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10971

GenBank Accession Number:

BC041622

GeneID (NCBI):

55734

UNIPROT ID:

Q9NTW7

Full Name:

zinc finger protein 64 homolog (mouse)

Calculated MW:

681 aa, 75 kDa

Purification Method:

Protein A purification

CloneNo.:

240847G6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00803-1, ZFP64 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83825-2-PBS. Detection antibody: 83825-1-PBS. Standard: Ag10971. Range: 0.313-40 ng/mL