

For Research Use Only

Mouse G-CSF Recombinant antibody, PBS Only (Detector)

Catalog Number: 83826-2-PBS



Basic Information

Catalog Number:

83826-2-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0256

GenBank Accession Number:

NM_009971

GeneID (NCBI):

12985

UNIPROT ID:

P09920

Full Name:

colony stimulating factor 3
(granulocyte)

Calculated MW:

22kd

Purification Method:

Protein A purification

CloneNo.:

240900D9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00806-1, MOUSE G-CSF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83826-1-PBS. Detection antibody: 83826-2-PBS. Standard: Eg0256. Range: 0.156-20 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83826-2-PBS against Mouse G-CSF was performed. The affinity constant is 2.19 nM.