

For Research Use Only

Mouse G-CSF Recombinant antibody, PBS Only (Detector)

Catalog Number:83826-4-PBS



Basic Information

Catalog Number: 83826-4-PBS	GenBank Accession Number: NM_009971	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 12985	CloneNo.: 250727G6
Source: Rabbit	UNIPROT ID: P09920	
Isotype: IgG	Full Name: colony stimulating factor 3 (granulocyte)	
Immunogen Catalog Number: EG0256	Calculated MW: 22kd	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00806-2, Mouse G-CSF Recombinant Matched Antibody Pair - PBS only. 83826-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0256.83826-4-PBS was HRP conjugated as the detection antibody. Range: 15.6-1000 pg/mL



Serum of sixteen mice was measured. The G-CSF concentration of detected samples was determined to be 594.3 pg/mL with a range of 189.4-987.2 pg/mL