

For Research Use Only

ATP2B1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83827-1-PBS



Basic Information

Catalog Number:

83827-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG32423

GenBank Accession Number:

NM_001001323.1

GeneID (NCBI):

490

UNIPROT ID:

P20020

Full Name:

ATPase, Ca⁺⁺ transporting, plasma membrane 1

Calculated MW:

135 kDa

Purification Method:

Protein A purification

CloneNo.:

240848C4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00808-2, ATP2B1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83827-1-PBS. Detection antibody: 83827-2-PBS. Standard: Ag32423. Range: 0.313-40 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83827-1-PBS against Human ATP2B1 was performed. The affinity constant is 22.3 nM.