

For Research Use Only

CLCA1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83829-3-PBS



Basic Information

Catalog Number:

83829-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18176

GenBank Accession Number:

BC156805

GeneID (NCBI):

1179

UNIPROT ID:

A8K7I4

Full Name:

chloride channel accessory 1

Calculated MW:

914 aa, 100 kDa

Purification Method:

Protein A purification

CloneNo.:

240941E6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00794-3, CLCA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83829-3-PBS. Detection antibody: 83829-6-PBS. Standard: Ag18176. Range: 0.625-80 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83829-3-PBS against Human CLCA1 was performed. The affinity constant is below 1 pM.