

For Research Use Only

# DLL3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83834-2-PBS



## Basic Information

Catalog Number:

83834-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0579

GenBank Accession Number:

BC000218

GeneID (NCBI):

10683

UNIPROT ID:

Q9NYJ7

Full Name:

delta-like 3 (Drosophila)

Calculated MW:

65 kDa

Purification Method:

Protein A purification

CloneNo.:

240956D5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00798-1, DLL3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83834-2-PBS. Detection antibody: 83834-1-PBS. Standard: Ag0579. Range: 0.313-40 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83834-2-PBS against Human DLL3 was performed. The affinity constant is 1.70 nM.