

For Research Use Only

# DLL3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83834-3-PBS



## Basic Information

Catalog Number:

83834-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0579

GenBank Accession Number:

BC000218

GeneID (NCBI):

10683

UNIPROT ID:

Q9NYJ7

Full Name:

delta-like 3 (Drosophila)

Calculated MW:

65 kDa

Purification Method:

Protein A purification

CloneNo.:

240956A11

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP00798-2, Human DLL3 Recombinant Matched Antibody Pair - PBS only. 83834-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag0579. 83834-4-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL

TT cell lysate was measured. The human DLL3 concentration of detected samples was determined to be 2.84 ng/mL (based on a 2.9 mg/mL extract load).