

For Research Use Only

# GOLGA5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83839-1-PBS

## Basic Information

Catalog Number:	83839-1-PBS	GenBank Accession Number:	BC023021	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	9950	CloneNo.:	240893E2
Source:	Rabbit	UNIPROT ID:	Q8TBA6		
Isotype:	IgG	Full Name:	golgi autoantigen, golgin subfamily a, 5		
Immunogen Catalog Number:	AG35031	Calculated MW:	83 kDa		

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00813-1, GOLGA5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83839-3-PBS. Detection antibody: 83839-1-PBS. Standard: Ag35031. Range: 0.313-40 ng/mL.