

For Research Use Only

GOLGA5 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83839-3-PBS

Basic Information

Catalog Number:

83839-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35031

GenBank Accession Number:

BC023021

GeneID (NCBI):

9950

UNIPROT ID:

Q8TBA6

Full Name:

golgi autoantigen, golgin subfamily a, 5

Calculated MW:

83 kDa

Purification Method:

Protein A purification

CloneNo.:

240893F7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00813-1, GOLGA5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83839-3-PBS. Detection antibody: 83839-1-PBS. Standard: Ag35031. Range: 0.313-40 ng/mL.