

For Research Use Only

GOLGA5 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83839-5-PBS

Basic Information

Catalog Number:	83839-5-PBS	GenBank Accession Number:	BC023021	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	9950	CloneNo.:	240893H8
Source:	Rabbit	UNIPROT ID:	Q8TBA6		
Isotype:	IgG	Full Name:	golgi autoantigen, golgin subfamily a, 5		
Immunogen Catalog Number:	AG35031	Calculated MW:	83 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00813-3, GOLGA5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83839-5-PBS. Detection antibody: 83839-2-PBS. Standard: Ag35031. Range: 0.313-40 ng/mL.