

For Research Use Only

# IQUB Recombinant antibody, PBS Only (Capture)



Catalog Number: 83847-4-PBS

## Basic Information

Catalog Number:	83847-4-PBS	GenBank Accession Number:	BC128181	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	154865	CloneNo.:	240957F10
Source:	Rabbit	UNIPROT ID:	Q8NA54		
Isotype:	IgG	Full Name:	IQ motif and ubiquitin domain containing		
Immunogen Catalog Number:	AG19507	Calculated MW:	791 aa, 93 kDa		

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00822-1, IQUB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83847-4-PBS. Detection antibody: 83847-1-PBS. Standard: Ag19507. Range: 0.625-80 ng/mL