

For Research Use Only

vwf Recombinant antibody, PBS Only (Capture)

Catalog Number: 83854-3-PBS



Basic Information

Catalog Number:

83854-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25578

GenBank Accession Number:

GeneID (NCBI):

7450

UNIPROT ID:

P04275

Full Name:

von Willebrand factor

Purification Method:

Protein A purification

CloneNo.:

240867H2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00810-1, VWF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83854-3-PBS. Detection antibody: 83854-1-PBS. Standard: Ag25578. Range: 0.313-40 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83854-3-PBS against Human vwf was performed. The affinity constant is 1.21 nM.