

For Research Use Only

BCHE Recombinant antibody, PBS Only (Capture)



Catalog Number: 83855-3-PBS

Basic Information

Catalog Number: 83855-3-PBS	GenBank Accession Number: BC008396	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 590	CloneNo.: 240863G6
Source: Rabbit	UNIPROT ID: P06276	
Isotype: IgG	Full Name: butyrylcholinesterase	
Immunogen Catalog Number: AG20888	Calculated MW: 602 aa, 68 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00804-2, BCHE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83855-3-PBS. Detection antibody: 83855-1-PBS. Standard: Ag20888. Range: 1.25-20 ng/mL.