

For Research Use Only

Mouse IL-9 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83877-2-PBS

Basic Information

Catalog Number:

83877-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_008373.1

GeneID (NCBI):

16198

UNIPROT ID:

P15247

Full Name:

interleukin 9

Calculated MW:

16 kDa

Purification Method:

Protein A purification

CloneNo.:

240818G1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00858-1, MOUSE IL-9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83877-2-PBS. Detection antibody: 83877-4-PBS. Standard: Eg0851. Range: 0.313-40 ng/mL