

For Research Use Only

# MRPL16 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83888-2-PBS

## Basic Information

Catalog Number: 83888-2-PBS	GenBank Accession Number: BC001040	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 54948	CloneNo.: 241069B2
Source: Rabbit	UNIPROT ID: Q9NX20	
Isotype: IgG	Full Name: mitochondrial ribosomal protein L16	
Immunogen Catalog Number: AG6523	Calculated MW: 28 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00832-2, MRPL16 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83888-2-PBS. Detection antibody: 83888-4-PBS. Standard: Ag6523. Range: 0.78-100 ng/mL