

For Research Use Only

MRPL16 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83888-4-PBS

Basic Information

Catalog Number: 83888-4-PBS	GenBank Accession Number: BC001040	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 54948	CloneNo.: 241069H8
Source: Rabbit	UNIPROT ID: Q9NX20	
Isotype: IgG	Full Name: mitochondrial ribosomal protein L16	
Immunogen Catalog Number: AG6523	Calculated MW: 28 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00832-2, MRPL16 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83888-2-PBS. Detection antibody: 83888-4-PBS. Standard: Ag6523. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00832-3, MRPL16 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83888-1-PBS. Detection antibody: 83888-4-PBS. Standard: Ag6523. Range: 0.78-100 ng/mL