

For Research Use Only

4EBP1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83890-3-PBS

Basic Information

Catalog Number:	83890-3-PBS	GenBank Accession Number:	BC004459	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	1978	CloneNo.:	241055G1
Source:	Rabbit	UNIPROT ID:	Q13541		
Isotype:	IgG	Full Name:	eukaryotic translation initiation factor 4E binding protein 1		
Immunogen Catalog Number:	AG5056	Calculated MW:	118 aa, 12 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00559-5, 4EBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83890-3-PBS. Detection antibody: 83890-2-PBS. Standard: Ag5056. Range: 0.625-80 ng/mL